




## Trident Western HRP Substrate Series

— Blot Better for Less —

Citation Support

Cat. No.	Product	package
GTX14698	Trident femto Western HRP Substrate	100/200µl
GTX17435	Trident pico Western HRP Substrate	100/500µl
GTX400006	Trident plus Western HRP Substrate	500µl

### Side by Side Comparison



WB analysis using Trident Western HRP Substrate Series with various sample amount and exposure times as indicated. SuperSignal™ West Femto, Clarity™, Luminata™ Forte are registered trademark of Thermo Fisher Scientific, Bio-Rad, and EMD Millipore, respectively. The trademark holder is not affiliated with GeneTex and does not endorse this product.

### Advantages



Excellent Stability



Strongest Signal



Lower Cost



Less Sample



Easy to use




## Protein Ladders

Citation Support

Cat. No.	Product	package
GTX50875	Trident Prestained Protein Ladder	500µl
GTX16376	Trident Blue Prestained Protein Ladder	500µl

### Trident Prestained Protein Ladder (GTX50875)

- ✓ 3-color protein standard with 12 pre-stained proteins
- ✓ 10 blue bands with 1 green (25 kDa) and 1 red (75 kDa) band
- ✓ Cover wide range molecular weights from 10 to 245 kDa
- ✓ Suitable and verified for Western transfer on PVDF, nylon, or nitrocellulose membranes




## Protein extraction kits

Cat. No.	Product	package
GTX400005	Trident RIPA Lysis Buffer	100 ml
GTX16372	Trident Total Protein Extraction Kit	5/20 Tests
GTX16373	Trident Membrane Protein Extraction Kit	5/20 Tests
GTX16374	Trident Nuclear Protein Extraction Kit	4/20 Tests

Trident Membrane protein extraction kit (**GTX16373**) is for rapid extraction of native total membrane proteins from cultured mammalian cells or tissues.

- ✓ Simple and user-friendly
- ✓ Wide range of starting cells (1 - 50 million / sample)
- ✓ Detergent and EDTA free
- ✓ No need for Dounce homogenizer or tissue blender
- ✓ Finish extraction in less than 45 minutes
- ✓ High yield



**SCAN FOR ONLINE PRODUCT DETAILS!**

